Integration, deconvolution or qQMSA ?

Reino Laatikainen and Bekker Laatikainen SpinDiscovenes Etd.

Examples by Hannu Maaheimo (V-FF) and Yujo Hiltunen (UEF)

Basic integration

When there is quantitative reference (TMS, TSP, DSS or REFi) and its concentration are given (in PMR-file), any spectral integrals (area) can be transformed into **mmol** when the number of protons in the integral area are defined. When the corresponding molecular weight (**MW**) is known, also **mg/ml** is obtained.







Example 2: kinetics by novel traditional integration

<u>Step 1</u>: normalize spectra so that the volume of **reference** (here TSP) is the same in all the spectra. <u>Step 2</u>: integrate **reference** signal and (if not TMS, TSP or DSS), name it, give no. of protons (9 inTSP) and molarity of it (MMOL). <u>Step 2</u>: select the signals to be integrated, name them, give no. of protons in each signal (for example, 3 for CH₃) and the MW of the corresponding compound. <u>Step 4</u>: Integrate all the spectra in stacked mode. <u>Step 5</u>: Show concentrations by ChemAdder tool, or transfer to EXCEL.

Easy, works well if signals well separated.

ChemAdder integration tool allows (spline polynomial) baseline fitting and therefore also small signals sitting on large broad signal can be safely integrated.

Integrate or fit (QMSA), then drag-and-drop to EXCEL or use the tools of ChemAdder – all data with same click!



Kinetics by qQMSA

<u>Step 1:</u> analyze one representetative spectrum to make QMSA model for the system, use ignore for extra the signals.

<u>Step 2:</u> batch analyze the spectra with QMSA.

<u>Step 3:</u> show the concentrations by ChemAdder tool or transfer to EXCEL.

Spectral parameters needed but works also if signals not separated.

Example 3: Analysis of a macromolecular material with impurities (example by Yrjö Hiltunen)



Deconvolution (TLS) analysis

Can be used to model signals which are not described by QM models – together with QMSA signals

The spectrum (not explained by QM model) can be composed from extra lines or regular multiplets (doublets, triplets, etc.), which can be also named and grouped. If one sets a signal range and the number of lines in it, SpinAdder creates a model assuming equidistant lines with same line-width – the positions and line-widths can then be optimized.



HOLISTIC qQMSA

<u>Combines qQMSA, TLS, integration, pattern search and forms mathematical</u> <u>model for correlations of metabolite chemical shifts</u>

See Presentation: ChemAdder_HOLISTICS

